REDLIGHT

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 **BLUETOUCH.**

**WARNING**

This device is not intended for diagnosis, treatment, or prevention of diseases. The company provides no explicit or implied guarantees regarding therapeutic effectiveness related to physical health, dermatological conditions, mental health, or any other medical field. Use of this product for therapeutic purposes must be supervised by qualified healthcare professionals and follow their instructions. Users should fully understand the intended purpose of this product. Any legal responsibility arising from misuse, health risks, or subsequent claims lies solely with the responsible party.

I. PRODUCT DESCRIPTION

The Red Light Therapy Cap is a device designed to support the treatment of hair loss and stimulate hair growth. It gently stimulates weakened hair follicles, increases their vitality, and promotes healthy hair growth. The device uses 850 nm infrared light and 660 nm red light.

II. INDICATIONS FOR USE

This product is intended for:

- supporting the treatment of thinning hair and baldness,
- helping slow down and prevent hair loss,
- stimulating hair regrowth.

III. CONTRAINDICATIONS

1. Do NOT use this device if you have any of the following:

- metastatic tumors,
- high fever,
- obliterative vasculitis,
- eye diseases,
- malignant tumors,
- vascular insufficiency,
- active tuberculosis,
- acute soft-tissue infections,
- pregnancy,
- bleeding disorders,
- skin sensitivity disorders,
- skin allergies.

2. Do NOT use if:

- you have sustained an acute injury within the last 72 hours,
- you have a pacemaker, defibrillator, or other electronic implant,
- you are photosensitive or taking medications that may cause photosensitivity. If any signs of photosensitivity occur, discontinue use immediately. Consult your physician if you take or may take medications that increase light sensitivity.

IV. WARNINGS

Do NOT use this device if you are:

- sleeping or unconscious,
- under the influence of alcohol or drugs.

Do NOT use the device if:

- you have reduced skin sensitivity due to medication or other causes,
- the cable or plug is damaged,
- the device is damaged, torn, or has come into contact with water,
- the device is near hot surfaces.

Do NOT:

- bend or fold the device,
- use it on open wounds (heat may increase bleeding),
- use it together with warming creams, gels, or topical agents (risk of skin burns),
- cover the device while it is operating.

This is an electrical device—handle with care.

V. PRECAUTIONS

Consult a doctor before use especially if you:

- are pregnant,
- have suspicious or cancerous skin lesions.

VI. PRODUCT DETAILS

The product is fully assembled and ready for use.

Package contents:

- Red Light Therapy Cap,
- USB power cable,
- User manual.

VII. LIGHT SENSITIVITY TEST

Before first use, perform the following test:

1. Turn on the device.
2. Place it on approximately 1/4 of your forearm for 3 minutes.
3. If strong redness appears and lasts longer than 2 hours, your skin is light-sensitive.
4. In this case, **do not use the device**.

VIII. INSTRUCTIONS FOR USE

- Press and hold the power button to turn the device on.
- Short presses change the light intensity: low / medium / high.
- The MODE button switches between modes:
 - continuous light,
 - flashing light,
 - 10 Hz flashing,
 - 40 Hz flashing.

Put the cap on your head.

Only the red LEDs will be visible—infrared LEDs emit light that is invisible to the human eye.

Recommended treatment time: 10–20 minutes daily.

To turn off the device, press and hold the power button.

The device automatically shuts off after 20 minutes.

IX. FUNCTION DESCRIPTION (BUTTONS SW1 AND SW2)

SW1 - Main Button

Long press: turn device on/off

- When turning on: 2 short beeps, red + IR lights turn on continuously.
- When turning off: 1 long beep.

Short presses:

- 1× - Low IR intensity,
- 2× - Medium IR intensity,
- 3× - High IR intensity.

SW2 - Mode Button

- 1× - Slow “breathing” effect (smooth pulsing),
- 2× - Fast flashing (10 Hz),
- 3× - Slower flashing (40 Hz),

- 4x – Continuous light.

Modes cycle repeatedly.

Charging Indicators:

- While charging: red light slowly flashes,
- Low battery: red light slowly flashes,
- Fully charged: green light stays on,
- After unplugging: indicator turns off.

The device stops functioning during charging.

X. CLEANING INSTRUCTIONS

Clean the LED surface once a month with isopropyl alcohol and a soft cotton cloth. Allow the device to cool before cleaning. After cleaning, wait until the alcohol has completely evaporated before using the device again.

XI. FREQUENTLY ASKED QUESTIONS (FAQ)

How often should I use the device?

Daily for 10–20 minutes for best results.

Do the LEDs need replacing?

No, they are designed to last for several years under normal use.

What light intensity should I choose?

We recommend always starting at the lowest intensity and increasing it according to your feelings and reactions.

XII. DISCLAIMER

This product is not intended to diagnose, treat, or prevent any disease. The manufacturer makes no medical claims regarding the effectiveness of this product. If you experience any health problems requiring medical attention, consult a qualified healthcare professional. The manufacturer is not responsible for safety issues resulting from improper use of the product.

XIII. TECHNICAL SPECIFICATIONS

All settings are factory-preset.

This product is not user-serviceable.

Operating temperature: 0 °C - 27 °C (32 °F - 80 °F)

Relative humidity: up to 85%

Storage temperature: -30 °C - +55 °C (-22 °F - +130 °F)

Wavelengths: IR 850 nm, Red 660 nm

Power supply: 5 V / 1 A

XIV. DEVICE CLASSIFICATION

Standard device (not protected against moisture ingress).

XV. LIMITED WARRANTY

The product must be used according to this manual and only with the supplied power cable.

XVI. AFTER-SALES SERVICE

Thank you for choosing our product.

If you encounter any issues during use, please contact our customer service team.

SAFETY INSTRUCTIONS

Safety Instructions

Use the device only for purposes described in the manual. Do not use if the device is visibly damaged, overheated, or emits an unusual smell. Persons with reduced heat sensitivity must use the device only under supervision. Do not use on open wounds, burned skin, fresh scars, infections, or inflamed areas. Do not use if you have an implanted electronic device (e.g., pacemaker) unless approved by a doctor. Keep out of reach of children unless stated otherwise. Do not use while sleeping. Do not expose the device to water or high humidity. Do not use near flammable materials.

Contraindications and Warnings

Not suitable for persons with photosensitive conditions. Not suitable for light-sensitive epilepsy. Not suitable during pregnancy (unless approved by a doctor). If irritation, redness, or pain occurs, stop using the device immediately.

Maintenance and Cleaning

Allow the device to cool down after use. Clean only with a soft dry or slightly damp cloth. Do not use alcohol, solvents, or aggressive cleaning agents. Do not immerse the device in water.

Storage

Store in a dry place away from direct sunlight. Storage temperature: -10 °C to +45 °C. Protect from dust and humidity.

Device Disposal

Do not dispose of with household waste. Take to a specialized collection point.

Manufacturer

BlueTouch, s.r.o.

Zelený pruh 1974/109,

140 00 Prague, Czechia

ID: 09437223, VAT: CZ09437223

info@bluetouch.cz

EU COMPLIANCE

This product has been designed and manufactured in accordance with the applicable legal regulations of the European Union relating to consumer electronics.

The product complies with the requirements of the following EU directives and regulations:

- Directive 2014/30/EU on electromagnetic compatibility (EMC),
- Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS),
- ČSN ISO 13485,
- ČSN EN IEC 55014-1.

Based on compliance with the above requirements, the product is marked with the CE marking.

The EU Declaration of Conformity is available from the importer and can be provided upon request to the relevant supervisory authorities or customers.

This product is intended exclusively for wellness and relaxation use.

It is not a medical device within the meaning of Regulation (EU) 2017/745 and is not intended for the diagnosis, treatment, or prevention of diseases.

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**MORE INFORMATION ABOUT LIGHT THERAPY
CAN BE FOUND ON OUR WEBSITE.**

BLUETOUCH.