REDLIGHT

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76 LED DIODES

OPERATING INSTRUCTIONS

1. Power on: Press and hold the power button for three seconds to turn on the device. If the device cannot be turned on normally, please charge it before using it. It is strictly forbidden to use it while charging.

2. Vibration settings: Tap the wave-shaped vibration level button in the lower left corner to switch vibration levels. Each tap switches the vibration level. There are three levels: M1 / M2 / M3.

3. Red light mode: Click the sun-shaped button in the middle to switch to red light mode. Each click switches the red light mode. There are five modes in total: L1 / L2 / L3 / L4 / L5.

4. Temperature settings: Click the temperature button in the lower right corner to switch the temperature. Each click will switch the temperature. There are three modes: 45 °C - 50 °C - 55 °C.

5. Lock (protection against accidental pressing): Press the power button once, a lock icon will appear on the screen, and the anti-accidental touch mode will be turned on. During this time, no buttons can be used.

6. Power off: When not in use, press and hold the power button for three seconds to turn off the device.

BASIC INFORMATION

Product name:	BLUETOUCH - light therapy for shoulder
LED diodes:	76 pcs
Voltage / Battery capacity:	5V 2A / 5 000 mAh
Wavelengths:	660 nm + 850 nm
Brightness levels:	L1 / L2 / L3 / L4 / L5
Vibration modes:	M1 / M2 / M3
Temperature:	45°C / 50°C / 55°C

WARNING

This device is not intended for diagnosis, treatment, or prevention of diseases. The company provides no explicit or implied guarantees regarding therapeutic effectiveness related to physical health, dermatological conditions, mental health, or any other medical field. Use of this product for therapeutic purposes must be supervised by qualified healthcare professionals and follow their instructions. Users should fully understand the intended purpose of this product. Any legal responsibility arising from misuse, health risks, or subsequent claims lies solely with the responsible party.

ADDITIONAL INFORMATION

ADDITIONAL INFORMATION

- The device automatically switches off after 30 minutes of use. To use it again, turn it back on.
- Do not use the device while it is charging.
- The device charges via a USB-C port. A charging cable is included.

SAFETY INSTRUCTIONS

Safety Instructions

Use the device only for purposes described in the manual. Do not use if the device is visibly damaged, overheated, or emits an unusual smell. Persons with reduced heat sensitivity must use the device only under supervision. Do not use on open wounds, burned skin, fresh scars, infections, or inflamed areas. Do not use if you have an implanted electronic device (e.g., pacemaker) unless approved by a doctor. Keep out of reach of children unless stated otherwise. Do not use while sleeping. Do not expose the device to water or high humidity. Do not use near flammable materials.

Contraindications and Warnings

Not suitable for persons with photosensitive conditions. Not suitable for light-sensitive epilepsy. Not suitable during pregnancy (unless approved by a doctor). If irritation, redness, or pain occurs, stop using the device immediately.

Maintenance and Cleaning

Allow the device to cool down after use. Clean only with a soft dry or slightly damp cloth. Do not use alcohol, solvents, or aggressive cleaning agents. Do not immerse the device in water.

Storage

Store in a dry place away from direct sunlight. Storage temperature: -10 °C to +45 °C. Protect from dust and humidity.

Device Disposal

Do not dispose of with household waste. Take to a specialized collection point.

Manufacturer

BlueTouch, s.r.o.
Zelený pruh 1974/109,
140 00 Prague, Czechia
ID: 09437223, VAT: CZ09437223
info@bluetouch.cz

EU COMPLIANCE

This product has been designed and manufactured in accordance with the applicable legal regulations of the European Union relating to consumer electronics.

The product complies with the requirements of the following EU directives and regulations:

- Directive 2014/30/EU on electromagnetic compatibility (EMC),
- Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS),
- ČSN ISO 13485,
- ČSN EN IEC 55014-1.

Based on compliance with the above requirements, the product is marked with the CE marking.

The EU Declaration of Conformity is available from the importer and can be provided upon request to the relevant supervisory authorities or customers.

This product is intended exclusively for wellness and relaxation use.

It is not a medical device within the meaning of Regulation (EU) 2017/745 and is not intended for the diagnosis, treatment, or prevention of diseases.

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**MORE INFORMATION ABOUT LIGHT THERAPY
CAN BE FOUND ON OUR WEBSITE.**

BLUETOUCH.